



# NEW PRODUCT GENERATION



HUR takes the next step in exercise intelligence launching already the 5th generation of computerized equipment. **HUR SmartTouch Hi5 strength training equipment** leads the way for smart exercise solutions in the fields of active ageing, senior exercise, rehabilitation and inclusive wellness. **The new equipment is designed to further enhance usability and training experience for selected user groups.**

- Larger screen, improving ease of use
- Android technology enabling new features and better user experience

# HUR SMARTTOUCH

## THE 5<sup>TH</sup> PRODUCT GENERATION



- Evolution of common technology platforms and tools
- Continuous product development and engineering
- Research work with science community

3rd Generation  
HUR SmartCard  
Year 2008



4th Generation  
HUR  
SmartTouch  
Year 2015



5th Generation  
HUR Hi5  
SmartTouch  
Year 2018



What's new...and latest technology in Hi5 SmartTouch?

- 10" Screen in portrait mode
- Advanced sensor technology
- Future platform – advanced features
- Electronics and mechanics re-engineered

# 4GEN VS. HI5 COMPARISON

- Wireless connection – no ethernet cables
- Bigger screen
- Platform for the upcoming years!
- New cushion type and one new color option
- New seat adjustment in Leg Extension/Curl



HUR SMARTTOUCH



# HUR HI5 MACHINES



HUR Senior/Rehab  
For Senior Exercise and Rehabilitation



9 machines \*

EASY ACCESS



9 machines \*



HUR Performance  
For Active silverfoxes,



12 Machines available 2018  
16 Machines available Q1 2019

HUR MULTIFUNCTION



2 Machines

PULLEY



4 Variants

\*Will be available as MD class I





# HUR PREMIUM LINE



## For rehabilitation and senior exercise

HUR's most developed range offers world-leading, university level equipment for use in senior exercise and rehabilitation. **In these sectors it is crucial, for efficient and safe training and operation, to have a computerized system for automated reporting and the possibility for individual safe training with pneumatic resistance** – as well as a close to zero starting load and 100g/1kg resistance increments.

This equipment range, which is suitable for all ages and abilities, incorporates many features that enhance safety and facilitate rehabilitation, for example, range limiters and additional support. All HUR Premium Line equipment is recognized as medical devices.

# HUR PREMIUM LINE



## Upper-body



5520-Hi5 Push up/ Pull down



5540-Hi5 Chest press



5160-Hi5 Peck deck



5175-Hi5 Optimal rhomb

## Core



5310-Hi5 Abdomen/ Back

## Easy Access (EA) Core



9330-Hi5 Twist EA



9310-Hi5 Abdomen, Back EA

## Lower-body



5520 -Hi5 Adduction/ Abduction



5510-Hi5 Body Extension



5530-Hi5 Leg Extension/ Curl



5540-Hi5 Leg Press

## Easy Access (EA) Upper-body



9175-Hi5 Optimal Rhomb EA



9160-Hi5 Peck deck EA



9150-Hi5 Lat Pull EA



9140-Hi5 Chest Press EA



9110-Hi5 Biceps/ Triceps EA



9125-Hi5 Dip/Shrug EA



9120-Hi5 Push up/ Pull down EA



## HUR ACTIVE LINE



### For active living

HUR's Active range offers easy to use machines with minimal adjustment points. The unique Natural Transmission system provides safe and effective workouts with **no damaging stress on joints and connecting tissues** – the workload is applied specifically to the muscles. Furthermore the silent machines take up very little space and add another design element to wellness centres, hotels, companies etc.

# HUR ACTIVE LINE



## Upper-body



3110-Hi5 Biceps/



3520-Hi5 Push up/ Pull down



3540-Hi5 Chest Press



3125-Hi5 Dip/ Shrug



3175-Hi5 Optimal Rhomb

## Core



3310-Hi5 Abdomen/Back Press



3320-Hi5 Back



3330-Hi5 Twist

## Lower-body



3540-Hi5 Leg Press



3530-Hi5 Leg Extension/Curl



3110-Hi5 Body Extension



3520-Hi5 Adduction/ Abduction

## Multifunction



8540-Hi5 Multifunction 4X



8530-Hi5 Multifunction 5X

## Pulleys



8820-Hi5 & 8830-Hi5 Wall Mounter



8821-Hi5 & 8831-Hi5 Free standing



8822-Hi5 & 8832-Hi5 2-1



8823-Hi5 & 8833-Hi5 3-1

### 2019 Machines

- Chestcurl
- Curlup
- Calfpress
- Incline
- Leg Press CC





# POWER & ROM



- Advanced features for Hi5 machines
- In Medica 2018 we will demonstrate Power per repetition and Range of Motion features.
- Power and ROM are commercially available during Q1 2019



# AVAILABILITY

- Deliveries will start from 1.1.2019
- Multifunction machines can be delivered in the beginning of February 2019


# MATERIALS

- Manual
- Installation instructions
- Product cards
- Price list
- Product catalog

# HUR SMARTTOUCH AND INTEGRATIONS

HUR SmartTouch can be integrated with other applications/programs using Web Services or by using an integration engine for the HL7 interface. **These solutions enable importing of clients/patients into HUR SmartTouch or exporting them to other applications.** Also the training information can be exported.

A local integration partner is needed to build the custom integration. For more information, contact HUR on [software@hur.fi](mailto:software@hur.fi).



## STE LEG/EXTENSION CURL - IMPROVEMENT

What's new?


- Back support is electrically adjusted (200mm range)

Availability

- Can be ordered after Medica 2018, will be delivered with first commercial Hi5 deliveries. Available for STE models of Leg/Extension Curl

Upgradability

- No possibility to upgrade old/existing machines, only for new production






## LEG PRESS INCLINE 5546

What's new?

- Redesigned lower frame for easier access
- 3545 and 5545 will be discontinued and replaced with 5546

Availability

- Estimated availability Q1/2019



	3545	5545	5546
Max Load	243 kg	153 kg	197 kg
Cylinders	2 x 63mm	2 x 50mm	1 x 80mm


## LEG PRESS CC - IMPROVEMENTS

What's new?

- Redesigned construction
- Adjustable start position

Availability

- Q2/2019



## 3300 TWIST - IMPROVEMENT

What's new?

- All models (9330 and 3330) have the same lower frame design

Availability

- Available after Medica



